

Zoom Stereo Microscope Model SZMN Series

Product Features :

- Head : Head inclined 45°, rotates 360°, binocular head or trinocular head.
- Inter-pupillary : Between 54-75mm.
- Eyepiece : WF10x/20
- Objective : Zoom objective magnification 0.7x – 4.5x
- Working distance : 100 mm
- Stand : B1 : Plain Focusing Stan, B2: Top & Bottom Halogen Stand, B3 : Top & Bottom LED Stand
- Option : Supp Lens 0.5x, 1.5x, 2x , Dioptor Eyepiece 15x, 20x, Retical Scale



SZMN-B1



SZMN-T1



SZMN-T3

● Accessory Diagram ●

