

AnyView Microscope

Adjustable focus and Magnification

2.0 M pixel
10X~200X

Leading Innovation>>>

*Experiencing Microcosm
Enjoying Vision Fun*

*Digital Microscope
The Microscope YOU NEED*




USB Digital Microscope

- Package contents:
- Anyview Digital Microscope
 - Device holder with 1/4" thread
 - 2 x transparent stand spacers
- System requirements:
- CPU Pentium/Celeron/AMD 1.0G or above
 - Windows XP SP2/Vista or 7
 - Available USB 2.0 or 1.1 port
 - VGA Card:16bit or above

Measurement functions.



Digital Microscope and Video Camera



Features:

- Magnification: 10x - 200x @ 17" monitor full screen view
- Field of View:
 - Min.: 1.9 mm x 2.5 mm @ 17 mm working distance (2.8 μm/pixel)
 - Max.: 100mm x 76 mm@300 mm working distance
- Image Sensor: 2.0 Mega pixels
- Video Capture Resolution 1600x1200, 1280x960, 640x480
- Color: 24bit RGB or YUV or MJPG
- Lens: High Precision Microscope Lens
- Focus Range: Manual from 8mm to 300mm
- White Balance: Auto
- Exposure time: Auto
- Light Source: Integrated adjustable LED illumination
- Length: 122 mm
- Diameter: 36 mm
- Output:
 - USB streaming video
- Power supply:
 - 5V DC from USB port

Application scope:

- Biological anatomy observation
- Industrial inspection and quality control
- Dermatology(skin/scalp) exam
- Printing and textile inspection
- School teaching and research tool
- Medical analysis tool
- Student trial tool
- Reading aid
- ...



Up



Be careful



Moisture