



Dual Stereo Microscope SFM 70 Series

SFM70 series is an newest dual Magnification Stereo microscope with changeable knob design. SFM 70series is an industrial grade, dual Magnification Stereo microscope. The Fluorescent mini ring light offers shadow-free illumination. Microscope head inclined 45, rotate 360, interpupillary distance 55-75 mm With eyepiece diopter adjustment. The SFM70 series is an economical choice where superior optical quality is a must, but zoom capability is not required.



SFM 70

Model	SFM7012	SFM7013	SFM7024
Eyepiece	WF10x	WF10x	WF10x
Objective	1x, 2x	1x, 3x	2x, 4x
Magnification	10x, 20x	10x, 30x	20x, 40x
Working distance (mm)	100	100	100