



Video Microscope ICS-305

Product Features

- Unbeatable and stereo image
- Easy image saving and processing
- Convenient Image- saving by PC, precise measurement and easy editing after inspection
- Contact type/Non-contact type
- Handy focus control
- Memory function in main unit
- Stage viewer that supports the most suitable illumination
- Stereo inspection at various angles
- Varies kind of Lens



ICS-P305B (Pen Type)

ICS-P305B (Gun Type)

Main Unit	Lamp		12V/100W Halogen Lamp
	Light Intensity		Adjustable (0-100W)
	Dimension		209(W) x 250(D) x 101(H)mm
	Power Supply		110/220 VAC, 50/60Hz
	Output Signal		Composite:1EA, S-Video:1EA
	Frame Capture	Memory Type	4 Frame
		Nomemory Type	-
Memory Display		4 Full Frame, quad display(Memory Type)	
Camera	Image Sensor		1/3" SONY Super HAD CCD™
	Pixel	NTSC	768(H)*494(V)
		PAL	752(H)*582(V)
	Horizontal Resolution		460 TV Line
	Signal Process		Digital Signal Process (DSP)
	AGC/White Balance		Auto
	Min Illumination		2 Lux
	ELC		ON
Power Cable & Optic Fiber		2M	
Stage	Illumination		3Way Side / Transmitted / Coaxial
	Focus		Fine / Coarse Focus
	Stage Moving		X:90mm * Y:70mm
	Dimension		400(W)*450(D)*335(H)mm
Lens	Magnification		X40,X80,X100,X200,X240,X300, X500,X600,X800,X1200,X2400